

May 17, 2024

## 2D differentiation of NPCs to neurons & astrocytes

This protocol is a draft, published without a DOI.

Jessie Buth<sup>1</sup>

<sup>1</sup>UCLA



Jessie Buth

UCLA

OPEN  ACCESS



**Protocol Citation:** Jessie Buth 2024. 2D differentiation of NPCs to neurons & astrocytes . [protocols.io https://protocols.io/view/2d-differentiation-of-npcs-to-neurons-amp-astrocyt-ddds226e](https://protocols.io/view/2d-differentiation-of-npcs-to-neurons-amp-astrocyt-ddds226e)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** May 08, 2024

**Last Modified:** May 17, 2024


**Protocol Integer ID:** 99474

### Abstract

Protocol to differentiate neural progenitor cells to mature neurons and astrocytes in 2D culture.

### Protocol materials

 Corning® Matrigel® Basement Membrane Matrix, LDEV-free, 10 mL **Corning Catalog #354234** Step 4

 Corning® Matrigel® Basement Membrane Matrix, Phenol Red-free, LDEV-free, 10 mL **Corning Catalog #356237**

Step 4



## Overview

### 1 Typical workflow:

d(-14) thaw the NPCs (in ROCKi)

- maintain in +EGF+FGF

d(-7) passage (in ROCKi)

- maintain in +EGF+FGF

d0 passage (in ROCKi)

d1 full media change to start differentiation to neurons/astrocytes

d7 passage (without ROCKi)

d14 passage (without ROCKi)

d42 Mitomycin C treatment (**for neurons only**)

~d60-100 replate neurons/astrocytes for experiments

#### Note

If more cells are needed, the NPCs can be maintained additional passages in +EGF+FGF before differentiation, but can't be refrozen and thawed again. If frozen/thawed again, they don't proliferate or differentiate well.

## Media Recipe

### 2 **NPC Base Media (500mL):**

240mL 50% DMEM/F12 (Hyclone VWR Cat # SH30023.02)

240mL 50% Neurobasal (FisherSci Cat # 21103-049)

2.5mL 0.5X N2 (100X; Invitrogen Cat # 17502048)

5mL 0.5X B27 without VA (50X; Invitrogen, 12587010)

5mL 1X Glutamax (100x Invitrogen Cat # 35050061)

5mL 1X NEAA (Invitrogen Cat # 11140050)

\*250 uL 5 ug/mL Human Insulin Solution (250mg per 500mL) (Sigma Cat # I9278-5ML)

1mL 1X Primocin (stock 500X, InvivoGen ant-pm-2)

50uL 1 µg/mL Heparin (10 mg/ml stock Sigma H3419-100KU)

**Note**

\*Look up the lot number of human insulin & make aliquots of 2.5mg.  
The concentration is typically 9-11 mg/mL -> 2.5 mg will be ~250uL.

**For NPC maintenance:**

+20 ng/mL bFGF (stock 100 ug/mL): 10uL/50mL NPC base

+20 ng/mL EGF (stock 200 ug/mL): 5uL/50mL NPC base

**For Neuron differentiation:**

d0-d80+

+20 ng/mL BDNF (stock 200 ug/mL): 5uL/50mL NPC base

+20 ng/mL GDNF (stock 100 ug/mL): 10uL/50mL NPC base

+1 uM cAMP (stock 10mM): 5uL/50mL NPC base

~d42 (can be d42-d49)

1hr pulse with 5 ug/mL Mitomycin C ready made solution (stock 10 mg/mL Sigma Cat# M5353)

^^ this is enrich for neurons, dividing cells die off or stop proliferating over the next 4-7 days

**For Astrocyte differentiation:**

d0-d14

+10 ng/mL hLIF (stock 100 ug/mL): 5uL/50mL

+10 ng/mL EGF (stock 200 ug/mL): 2.5uL/50mL

d14-d80+

+10 ng/mL hLIF (stock 100 ug/mL): 5uL/50mL

+10 ng/mL CNTF (stock 100 ug/mL): 5uL/50mL

## Thawing Neural Progenitor Cells (NPCs)

- 3 See protocol for thawing hematopoietic progenitor cells (HPCs).  
<https://www.protocols.io/view/thawing-frozen-hematopoietic-stem-cells-hpcs-dddt226n>

**It is the same procedure for NPCs, except:**

- the plates are coated with **0.5 mg/mL standard matrigel**
- the cells are maintained in **NPC base +EGF+FGF+10uM ROCKi** (ROCKi is removed 24hr after thawing)
- the cells are plated at **250-300K/well in a 6w plate**

## Passaging NPCs (every 6-7 days - volumes listed at for 1w of a 6w plate)





## 4 Before starting, coat 6-well plates with 0.5 mg/mL standard matrigel (Set aside for later)

30m

## Note

**Use standard matrigel (regular or phenol-free):**









 Corning® Matrigel® Basement Membrane Matrix, LDEV-free, 10 mL **Corning Catalog #354234**

 Corning® Matrigel® Basement Membrane Matrix, Phenol Red-free, LDEV-free, 10 mL **Corning Catalog #356237**

**To make aliquots of matrigel:**

- The concentration of a lot of matrigel can be looked up here in "quality certificate lookup" section. (<https://www.corning.com/worldwide/en/products/life-sciences/resource-library.html>).
- To keep coating consistent across lots, I check the concentration and make aliquots. However, if you don't have the lot number information, 0.5mg/mL is close to a 1:20 dilution.

**To coat plates:**






- Thaw an aliquot of matrigel at  4 °C or  On ice
- Dilute with the appropriate amount of ice cold DMEM/F12  
 6.5 mL for 6.5-wells = 3.25 mg matrigel
- Aliquot  1 mL per well in a cold 6-well plate
- Shake/tilt the plate until the matrigel evenly coats the bottom
- Incubate  37 °C for  00:30:00
- After incubating, aspirate the matrigel
- Add  2 mL per well DMEM/F12
- Place the plate in the incubator at  37 °C

## Note

Plates can be stored in the incubator with DMEM/F12 for up to 1 week before use

4.1 Pre-treat the cells with  10 micromolar (μM) ROCK inhibitor (ROCKi)

1h 30m


- Pick up the media in the well to a 15mL conical tube.
- Add  10 micromolar (μM) ROCKi (  2 μL per well ) and mix well.
- Gently return the media to the dish.
- Incubate the cells with ROCKi for  00:30:00 -  01:00:00 at  37 °C



#### 4.2 Lift the cells with Accutase.

- Aspirate the media.
- Wash with 2mL/well PBS (without Ca/Mg).
- Aspirate the PBS.
- Add 1mL/well Accutase.
- Incubate 5 minutes at 37C.
- Pick up the cells in accutase with a glass pipette and blow the cells off the dish.
- Transfer the cells to a 15mL conical tube.
- Add 6mL NPC base media to the dish and blow the remaining cells off the dish.
- Transfer the remaining cells to the 15mL conical tube.
- Pipette up/down 3-5x in the conical tube to break apart clumps. To break apart clumps, pipette the cells up, and then press the pipette against the side of the tube and blow out quickly.
- Pick up the entire cell mixture and filter through a 100um cell strainer over a 50mL conical tube to remove any remaining large clumps.
- Transfer the cell mixture back to the 15mL conical tube.
- Centrifuge at 180g for 5 minutes room temperature.
- Aspirate the media.
- Close the conical tube lid and tap the pellet to loosen the cell pellet.
- Resuspend in 2-3mL of NPC base media.
- Pipette up/down 3-5x or until the cell solution is cloudy.

#### 4.3 Count the cells and aliquot to new plates.

- Count 10uL cells + 10uL trypan blue with a hemocytometer.
- Determine the #uL for 1 million cells.
- Add 1 million cells/well in  2 mL per well NPC+EGF+FGF+ROCKi in new 0.5 mg/mL matrigel coated plates.

**NPC+EGF+FGF+ROCKI (10mL):**

10mL NPC base media

1uL 20ng/mL EGF (stock 200ug/mL)

2uL 20ng/mL FGF (stock 100ug/mL)

10uL 10uM ROCKi (stock 10mM)

- Gently place the plate in the incubator and shake evenly (4x forward/back, 4x side-to-side, 4x forward/back).

#### Note

If the cells were confluent, 1 million cells/well is typically around a 1:2-1:4 passage dilution.



## Expansion of NPCs

- 5 d1 after passaging NPCs (see STEP 3)
- Perform a full media change and replace with NPC maintenance media without ROCKi (NPC base +EGF+FGF)  
**NPC+EGF+FGF (10mL):**  
10mL NPC base media  
1uL 20ng/mL EGF (stock 200ug/mL)  
2uL 20ng/mL FGF (stock 100ug/mL)
- d3
- Perform a half media change. Remove 900uL/well, and add a fresh 1mL/well NPC base +EGF+FGF.
- d5
- Perform a half media change. Remove 900uL/well, and add a fresh 1mL/well NPC base +EGF+FGF.
- d6-7
- Passage to new dishes coated with 0.5mg/mL standard matrigel in NPC base +EGF+FGF.  
for 6w plates -> 1 million cells/well in 2mL NPC base +EGF+FGF  
for 5cm cellBIND dishes -> 2.4 million cells/dish in 6mL NPC base +EGF+FGF

### Note

The high plating density is important, the cells need to be very dense to maintain progenitor identity.

I've stained the cells after maintaining in +EGF+FGF at high density up to 5x passages after thawing and > 90% of the cells remain PAX6+SOX2+NESTIN+. If passaged at lower density, the cells will stop dividing and/or differentiate.

## Differentiation of NPCs to Neurons

- 6 d1 after passaging NPCs (see STEP 3)
- Perform a full media change and replace with neuron differentiation media (NPC base +BDNF+GDNF+cAMP)  
**NPC+BDNF+GDNF+cAMP (10mL):**  
10mL NPC base media  
1uL 20 ng/mL BDNF (stock 200 ug/mL)  
2uL 20 ng/mL GDNF (stock 100 ug/mL)  
1uL 1 uM cAMP (stock 10mM)
- d3



- Perform a half media change. Remove 900uL/well, add a fresh 1mL/well NPC base +BDNF+GDNF+cAMP.

d5

- Perform a half media change. Remove 900uL/well, add a fresh 1mL/well NPC base +BDNF+GDNF+cAMP.

d7 [*differentiation passage 1 - no ROCKi during this step*]

- Passage to new dishes coated with 0.5mg/mL standard matrigel in NPC base +BDNF+GDNF+cAMP.
  - for 6w plates -> 1 million cells/well in 2mL media
  - for 5cm cellBIND -> 2.4 million cells/dish in 6mL media
  - for 10cm cellBIND dishes -> 5.8 million cells/dish in 12mL media

d9

- Perform a full media change and replace with neuron differentiation media (NPC base +BDNF+GDNF+cAMP) to remove any dead cells.

d11

- Perform a half media change. Remove 900uL/well, add a fresh 1mL/well NPC base +BDNF+GDNF+cAMP.

d13 [*differentiation passage 2 - no ROCKi during this step*]

- Passage to new dishes coated with 0.5mg/mL standard matrigel in NPC base +BDNF+GDNF+cAMP.
  - for 6w plates -> 400K cells/well in 2mL media
  - for 5cm cellBIND -> 1 million cells/dish in 6mL media
  - for 10cm cellBIND dishes -> 2.5 million cells/dish in 12mL media

d16+

- Perform a half media change every 3-4 days until 90-95% confluent. When the cells are 90-95% confluent, passage at the same density as d13 and continue with half media changes every 3-4 days until 90-95% confluent again.
  - for 5cm cellBIND -> remove ~2-2.5mL (depending how much media evaporated, leave enough media to cover the bottom), and add a fresh 3mL/well NPC base +BDNF+GDNF+cAMP.
  - for 10cm cellBIND -> remove ~4.5-5mL (depending how much media evaporated, leave enough media to cover the bottom), and add a fresh 6mL/well NPC base +BDNF+GDNF+cAMP.

**Note**

Seeding at the lower density and feeding the cells less often (every 4 days) results in more progenitor differentiation, plan to have around the number of dishes you want to end with after the 2nd differentiation passage.

~d42

- To enrich for neurons, around d42 (or close to this date when the dishes are 95% confluent again), treat the cells with 5ug/mL Mitomycin C for 1hour at 37C.
- After 1hr, perform 3x washes with 6mL PBS, and then replace with 6mL (5cm dishes) or 12mL (10cm dish) of fresh NPC base +BDNF+GDNF+cAMP.

**Note**

If the cells were ~95% confluent before Mitomycin C treatment, 1x 10cm typically yields around 3-6 million cells.

Continue performing half media changes until used for experiments. After Mitomycin C treatment, the cells should not expand further, so they likely won't need to be passaged again until plating for experiments.

*I've been using these for experiments ~d60-100.*

## Differentiation of NPCs to Astrocytes

- 7 d1 after passaging NPCs (see STEP 3)
- Perform a full media change and replace with astrocyte differentiation media #1 (NPC base +hLIF+EGF)  
NPC+hLIF+EGF (10mL):  
10mL NPC base media  
1uL 10 ng/mL hLIF (stock 100 ug/mL)  
0.5uL 10 ng/mL EGF (stock 200 ug/mL)
- d3
- Perform a half media change. Remove 900uL/well, and add a fresh 1mL/well NPC base +hLIF+EGF.
- d5
- Perform a half media change. Remove 900uL/well, and add a fresh 1mL/well NPC base +hLIF+EGF.
- d6-7 [*differentiation passage 1 - no ROCKi during this step*]
- Passage to new dishes coated with 0.5mg/mL standard matrigel in NPC base +hLIF+EGF.





for 6w plates -> 1 million cells/well in 2mL media  
for 5cm cellBIND -> 2.4 million cells/dish in 6mL media  
for 10cm cellBIND dishes -> 5.8 million cells/dish in 12mL media

~d9

- Perform a full media change and replace with neuron differentiation media (NPC base +hLIF+EGF) to remove any dead cells.

~d11

- Perform a half media change. Remove 900uL/well, and add a fresh 1mL/well NPC base +hLIF+EGF.

~d13-14 [*differentiation passage 2 - no ROCKi during this step*]

- Passage to new dishes coated with 0.5mg/mL standard matrigel in **NPC base +hLIF+CNTF**.  
for 6w plates -> 400K cells/well in 2mL media  
for 5cm cellBIND -> 1 million cells/dish in 6mL media  
for 10cm cellBIND dishes -> 2.5 million cells/dish in 12mL media

**NPC+hLIF+CNTF (10mL):**

10mL NPC base media

1uL 10 ng/mL hLIF (stock 100 ug/mL)

1uL 10 ng/mL CNTF (stock 100 ug/mL)

~d16+

- Perform a half media change every 3-4 days with **NPC base +hLIF+CNTF** until 95-100% confluent.
- When the cells are 95-100% confluent, passage at the same density as ~d13 and continue with half media changes with **NPC base +hLIF+CNTF** every 3-4 days until 95-100% confluent again.

for 5cm cellBIND -> remove ~2-2.5mL (depending how much media evaporated, leave enough media

to cover the bottom), and add a fresh 3mL/well **NPC base +hLIF+CNTF**.

for 10cm cellBIND -> remove ~4.5-5mL (depending how much media evaporated, leave enough media

to cover the bottom), and add a fresh 6mL/well **NPC base +hLIF+CNTF**.



### Note

If the cells were ~95% confluent before passaging, 1x 10cm typically yields around 9-14 million cells.

The cells will expand slower with each passage, avoid passaging until 95-100% confluent. Continue performing half media changes until used for experiments, around d50-60 astrocytes with more complex branching should appear and increase in number.

*I've been using these for experiments ~d60-100.*

*After replating for experiments or if maintained past d100, I change the media to NPC +BDNF+GDNF+cAMP to prevent the astrocytes from overgrowing and detaching from the plate.*